Exercise 6

Sugar analysis in honey

Example of calculations

A. <u>Direct analysis of the reducing sugars (glucose</u> & fructose) in honey-Layne-Eynon method

Process

- 1. Reduction of Fehling's solution (Fel. A + Fel. B) with standard invert sugar (glucose + fructose 1:1) solution
 - Determination of Fehlings solution "title" (i.e. the amount of invert sugar in grams required to reduce the amount of Cu⁺² contained in the titrated Fehling's solution)
- 2. Reduction of Fehling's solution with the honey solution (1 %w/v honey in water)
 - Determination of <u>the reducing sugars</u> in honey (glucose and fructose)

A. Direct analysis of the reducing sugars (glucose & fructose) in honey-Layne-Eynon method

Process

1. Test titrations result: e.g. a ml of standard or honey solution are consumed


2. Final titrations result: e.g. b ml of standard or honey solution are consumed

3. Calculations are done using the result of the final titrations (b ml)

Principle

- Hydrolysis of the sucrose (which is a non-reducing sugar) contained in the sample, at 60 °C, with HCl sol., to the reducing sugars glucose & fructose
- Determination of the total reducing sugars
 by the Layne-Eynon method as previously
- Calculation of sucrose by the equation:

(total reducing sugars - directly assayed reducing sugars) \times 0,95



Glucose + fructose

sucrose

Calculations example

- Standard solution of invert sugar: initial concentration 1% w/v invert sugar → diluted as: 50 ml in 200 ml with water → therefore the final concentration of the standard solution is 0,25 %w/v invert sugar
 - ❖ for the Layne-Eynon method, 42 ml of this solution were consumed to reduce the Fehling's solution
- Honey sample for direct analysis of reducing sugars (Honey 1): prepare a solution of 1% w/v honey in water
 - ❖ for the Layne-Eynon method, 13 ml of this solution were consumed to reduce the Fehling's solution
- Honey sample for analysis of total sugars (Honey 2): dilute 50 ml of Honey 1 to 100 ml and hydrolyze the contained sugar by HCl → the final concentration of this solution is 0,5 %w/v honey in water
 - ❖ for the Layne-Eynon method, 24 ml of this solution were consumed to reduce the Fehling's solution

Calculations example

- Calculation of the Fehling's solution "title":
 - ▶ 100 ml of standard solution contain 0,25 g of reducing sugars (invert sugar) the 42 ml that were consumed X g?? \rightarrow X = $(42 \times 0,25) / 100 = <math>0,105$ g reducing sugars (so, 0.105 g invert sugar is the "title" of the Fehling's solution!)
- Direct determination of reducing sugars in hioney:
 - ► The 13 ml of Honey 1 (1%w/v) that were consumed to reduce the Fehling's solution must contain 0,105 g reducing sugars

```
100 ml of Honey 1 sol. X?? \rightarrow X = (100 × 1,105) / 13 = 0,808 g of reducing sugars
```

- ▶ But the 100 ml of Honey 1 solution contain 1 g of honey Therefore, 0,808 g of reducing sugars are contained in 1 g of honey so, in 100 g of honey X?? \rightarrow X = (100 × 0.808) / 1 = 80,8 g of reducing sugar
- ► Therefore the honey contains 80,8 % w/w reducing sugars

Calculations example

- Determination of sucrose in honey:
 - The 24 ml of Honey 2 solution (0,5%w/v honey in water) that were consumed, contain 0,105 g of invert sugar
 The 100 ml of Honey 2 sol. X?? → X = (100 × 1,105) / 24 = 0,438 g of total reducing sugars (since any contained sucrose was hydrolyzed to glucose +fructose)
 - ▶ But 100 ml of of Honey 2 solution contain 0,5 g of honey therefore, in 0,5 g of honey, 0,438 g of total reducing sugars are contained In 100 g of honey X?? → X = (100 × 0.438) / 0,5 = 87,5 g of total reducing sugars
- ► Therefore honey contains: $(87,5-80,8) \times 0,95 = 6,37 \%$ w/w sucrose